


Product Datasheet

VPS45 antibody (orb19289)


Description	Goat polyclonal antibody to VPS45
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	VPS45 (C Terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 60kDa band observed in Human Testis and Peripheral Blood Mononucleocyte lysates (calculated MW of 65.1kDa according to NP_009190.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	ESSQVTSRSASRR
Clonality	Polyclonal
MW	65.1
Entrez	11311
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb19289 (0.3ug/ml)
staining of Human
Pe...



Immunohistochemical
staining of Human
Co...



Immunohistochemical
staining of Human
Pl...