



Product Datasheet

STXBP1 antibody (orb19291)

Description	Goat polyclonal antibody to STXBP1
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	MUNC18 / STXBP1 (isoform a)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 75kDa band observed in Human and Mouse Brain lysates (calculated MW of 68.7kDa according to NP_003156.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	SRVSFEDQAPTME
Clonality	Polyclonal
MW	68.7
Entrez	6812
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Brain lys...



Immunohistochemical staining of Human Co...



Western blot analysis of Human Cerebellu...