


Product Datasheet


FXR1 antibody (orb19295)

Description	Goat polyclonal antibody to FXR1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	FXR1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 80kDa band observed in lysates of cell line NIH/3T3 (calculated MW of 59.9kDa according to NP_001013457.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Khandjian et al Hum Mol Genet. 1998 Dec;7(13):2121-8.; PMID: 9817930). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	RIEGDNENKLPRED
Clonality	Polyclonal
MW	69.7; 60.8; 59.9
Entrez	8087
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.




orb19295 (0.2ug/ml)
staining of NIH/3T3

...



WB analysis of NIH-3T3 cell lysate using...



Immunofluorescence analysis of A431 cell...