


Product Datasheet

DSCAM antibody (orb19296)

Description	Goat polyclonal antibody to DSCAM
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	DSCAM
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 170kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 174kDa according to NP_996770.1). Recommended concentration: 2-6 µg/ml.
Protein Sequence	DSWDSAQRTKDVSPQ
Clonality	Polyclonal
MW	222; 220
Entrez	1826
Dilution Range	ELISA: 1:16000, WB: 2-3 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of paraforma...



Western blot analysis of HepG2 cell lysa...