

Product Datasheet

IFNGR2 Antibody (C-term) (orb1929666)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	This IFNGR2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 308-337 amino acids from the C-terminal region of human IFNGR2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37806 Da
Uniprot ID	P38484
NCBI	NP_005525.2
Dilution Range	WB: 1:1000, WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.