

Product Datasheet

HOMER3 Antibody (Center) (orb1929669)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	This HOMER3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 281-308 amino acids from the Central region of human HOMER3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39836 Da
Uniprot ID	Q9NSC5
NCBI	NP_004829.3, NP_001139194.1, NP_001139193.1, NP_001139196.1
Dilution Range	WB: 1:1000, WB: 1:1000, WB: 1:1000, IHC-P: 1:10~50, FC: 1:10~50
Expiration Date	12 months from date of receipt.