

Product Datasheet

SELL Antibody (C-term) (orb1929671)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	This SELL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 346-372 amino acids from the C-terminal region of human SELL.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42187 Da
Uniprot ID	P14151
NCBI	NP_000646.2
Dilution Range	WB: 1:2000, IHC-P: 1:10~50, FC: 1:25
Expiration Date	12 months from date of receipt.