

Product Datasheet

LAG3 Antibody (Center) (orb1929672)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	This LAG3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 103-132 amino acids from the Central region of human LAG3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57449 Da
Uniprot ID	P18627
NCBI	NP_002277.4
Dilution Range	WB: 1:2000, WB: 1:2000, WB: 1:2000, WB: 1:2000, WB: 1:1000, WB: 1:2000, WB: 1:2000
Expiration Date	12 months from date of receipt.