

Product Datasheet

HMHA1 Antibody (Center) (orb1929674)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	This HMHA1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 877-903 amino acids from the Central region of human HMHA1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	124614 Da
Uniprot ID	Q92619
NCBI	NP_001245257.1, NP_001269264.1, NP_001269263.1, NP_036424.2
Dilution Range	WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.