


Product Datasheet

CHRM1 antibody (orb19298)

Description	Goat polyclonal antibody to CHRM1
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	M1 mAChR / CHRM1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 48kDa band observed in Human Brain lysates (calculated MW of 51.4kDa according to NP_000729.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	DKRRWRKIPKRPGSVH
Clonality	Polyclonal
MW	51.4
Entrez	1128
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Brain lys...



Western blot analysis of Mouse Heart lys...