

Product Datasheet

DOM3Z Antibody (N-term Y88) (orb1929925)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	This DOM3Z antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 73-99 amino acids from the N-terminal region of human DOM3Z.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	44929 Da
Uniprot ID	O77932
NCBI	NP_005501.2
Dilution Range	WB: 1:1000, WB: 1:1000, FC: 1:10~50
Expiration Date	12 months from date of receipt.