

Product Datasheet

Calcipressin 3 (DSCR1L2) Antibody (N-term) (orb1930096)

Description

Purified Rabbit Polyclonal Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications
Target

This Calcipressin 3 (DSCR1L2) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 19-51 amino acids from the N-terminal region of human Calcipressin 3 (DSCR1L2).

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

27492 Da

Uniprot ID
Q9UKA8
NCBI

NP_001238911.1, NP_001238912.1, NP_038469.1, NP_001238908.1, NP_001238909.1, NP_001238907.1, NP_001238910.1, NP_001238914.1, NP_001238906.1

Dilution Range

WB: 1:1000

Expiration Date

12 months from date of receipt.