

Product Datasheet

M Sirt3 Antibody (C-term) (orb1930161)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This Sirt3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 304-334 amino acids from the C-terminal region of Mouse Sirt3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36615 Da
Uniprot ID	Q8R104
NCBI	NP_001171275.1 , NP_071878.2 , NP_001120823.1
Dilution Range	WB: 1:2000, WB: 1:64000, WB: 1:2000, WB: 1:2000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.