

Product Datasheet

RAI17 Antibody (C-term) (orb1930164)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This RAI17 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1038-1067 amino acids from the C-terminal region of human RAI17.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	115483 Da
Uniprot ID	Q9ULJ6
NCBI	NP_065071.1
Dilution Range	WB: 1:2000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.