

## Product Datasheet

### SMYD5 (RAI15) Antibody (C-term) (orb1930165)

**Description** Purified Rabbit Polyclonal Antibody (Pab)

**Species/Host** Rabbit

**Reactivity** Human, Mouse

**Conjugation** Unconjugated

**Tested Applications** IHC-P, WB

**Target** This SMYD5 (RAI15) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 341-371 amino acids from the C-terminal region of human SMYD5 (RAI15).

**Form/Appearance** Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.

**Note** For research use only

**Isotype** Rabbit IgG

**Clonality** Polyclonal

**MW** 47341 Da

**Uniprot ID** [Q6GMV2](#)

**NCBI** [NP\\_006053.2](#)

**Dilution Range** WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100

**Expiration Date** 12 months from date of receipt.