

Product Datasheet

OASL Antibody (C-term) (orb1930167)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This OASL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 484-514 amino acids from the C-terminal region of human OASL.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59226 Da
Uniprot ID	Q15646
NCBI	NP_937856.1 , NP_003724.1 , NP_001248754.1
Dilution Range	WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.