

Product Datasheet

OASIS Antibody (C-term) (orb1930168)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This OASIS antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 491-519 amino acids from the C-terminal region of human OASIS.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57005 Da
Uniprot ID	Q96BA8
NCBI	NP_443086.1
Dilution Range	WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.