

Product Datasheet

STEAP2 Antibody (Center) (orb1930441)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	This STEAP2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 229-258 amino acids of human STEAP2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56056 Da
Uniprot ID	Q8NFT2
NCBI	NP_001035756.1, NP_001231875.1, NP_001231873.1, NP_694544.2, NP_001035755.1
Dilution Range	WB: 1:500-1:1000, WB: 1:500, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.