

Product Datasheet

PAK1 antibody (orb19307)

Description	Goat polyclonal antibody to PAK1
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PAK1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 70kDa band observed in Human, Mouse and Rat Brain lysates (calculated MW of 60.6kDa according to NP_002567.3). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NTEKQKKKPKMSDE
Clonality	Polyclonal
MW	60.6
Entrez	5058
Dilution Range	ELISA: 1:1000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Brain lysat...