

Product Datasheet

SC24A Antibody (Center) (orb1930853)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This SC24A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 636-663 amino acids from the Central region of human SC24A.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	119749 Da
Uniprot ID	O95486
NCBI	NP_068817.1, NP_001239160.1
Dilution Range	WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.