

Product Datasheet

QKI antibody (orb19312)

Description	Goat polyclonal antibody to QKI
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Quaking / QKI
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 38kDa band observed in mouse brain lysates (calculated MW of 35.8kDa according to NP_068681.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	ETKEKPKTPDY
Clonality	Polyclonal
MW	37.7; 35.1; 36.0; 35.2
Entrez	9444
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19312
(0.3ug/ml)
staining of
mouse br...