

Product Datasheet

Phospho-ERK3-S189 Antibody (orb1931271)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Target	This Phospho-ERK3-S189 antibody is generated from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding S189 of human ERK3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	82681 Da
Uniprot ID	Q16659
NCBI	NP_002739.1
Dilution Range	IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.