



Product Datasheet

VHL antibody (orb19313)

Description	Goat polyclonal antibody to VHL
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	VHL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 17kDa band observed in Human Ovary and Testis lysates and in lysates of all cell lines tested (calculated MW of 19.7kDa according to NP_937799.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	RSLVKPENYRRLD
Clonality	Polyclonal
MW	24.2; 19.7
Entrez	7428
Dilution Range	ELISA: 1:128000, WB: 0.01-0.03 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Ovary lys...



Immunohistochemical staining of Human Ki...