

Product Datasheet

LGR5 (GPR49) Antibody (Center) (orb1931495)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	This LGR5 (GPR49) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 443-473 amino acids from the Central region of human LGR5 (GPR49).
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	99998 Da
Uniprot ID	O75473
NCBI	NP_001264155.1 , NP_003658.1
Dilution Range	WB: 1:1000, WB: 1:1000, IHC: 1:100, IHC-P: 1:10~50, FC: 1:10~50
Expiration Date	12 months from date of receipt.