

Product Datasheet

SULT2A Antibody (C-term) (orb1931560)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This SULT2A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 253-285 amino acids from the C-terminal region of human SULT2A.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33780 Da
Uniprot ID	Q06520
NCBI	NP_003158.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.