

Product Datasheet

FBLN5 antibody (orb19318)

Description	Goat polyclonal antibody to FBLN5
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Fibulin 5 / FBLN5
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 50kDa band observed in Human Colon, Heart and Ovary lysates (calculated MW of 50.2kDa according to NP_006320.2). Recommended concentration: 0.1-0.3 µg/ml. An additional band of unknown identity was also consistently observed at 23kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported
Protein Sequence	RPIKGPREIQLDLE
Clonality	Polyclonal
MW	50.2
Entrez	10516
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19318
(0.1ug/ml)
staining of
Human Ov...