

Product Datasheet

HIST1H2BJ/HIST1H2BK/HIST3H2BB Antibody(C-term) (orb1932102)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This HIST1H2BJ/HIST1H2BK/HIST3H2BB antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 98-126 amino acids from the C-terminal region of human HIST1H2BJ/HIST1H2BK/HIST3H2BB.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	13904 Da
Uniprot ID	P06899
NCBI	NP_066402.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.