

Product Datasheet

LGALS3BP Antibody (Center) (orb1932111)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This LGALS3BP antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 383-412 amino acids from the Central region of human LGALS3BP.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65331 Da
Uniprot ID	Q08380
NCBI	NP_005558.1
Dilution Range	WB: 1:2000, WB: 1:2000, IHC-P-Leica: 1:100
Expiration Date	12 months from date of receipt.