


Product Datasheet

CCM2 antibody (orb19322)

Description	Goat polyclonal antibody to CCM2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	CCM2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 48kDa band observed in Human Heart lysates (calculated MW of 48.8kDa according to NP_113631.1). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	KGEKSRDKKAHEK
Clonality	Polyclonal
MW	51.4; 48.8; 39.5
Entrez	83605
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



orb19322
(0.03ug/ml)
staining of
Human H...