



Product Datasheet

FBXL3 antibody (orb19323)

Description	Goat polyclonal antibody to FBXL3
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	FBXL3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 50-55kDa band observed in Human Adipose lysates (calculated MW of 48.7kDa according to NP_036290.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	KRGGRDSDRNSSEE
Clonality	Polyclonal
MW	48.7
Entrez	26224
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Adipose L...



Immunohistochemical staining of Human Th...