



Product Datasheet

GFAP antibody (orb19324)

Description	Goat polyclonal to GFAP
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GFAP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.IHC: In paraffin embedded Human Cerebellum shows specific staining of astrocytes. Recommended concentration, 2-4µg/ml.WB: Sandwich-type ELISA with increasing amount of recombinant GFAP captured by a rabbit antibody. Recommended reporter concentration: 5-10ug/ml. Please note this data was obtained using a previous batch.
Protein Sequence	DGEVIKESKQEHKD
Clonality	Polyclonal
Entrez	2670
Dilution Range	ELISA: 1:128000, WB: 0.001-0.1 µg/ml
Expiration Date	12 months from date of receipt.



orb19324
(0.001µg/ml)
staining of
Human...



Western blot
analysis of
Human Brain
(Ce...