

Product Datasheet

SNX26 antibody (orb19325)

Description

Goat polyclonal antibody to SNX26


Species/Host

Goat

Reactivity


Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications


orb19325
1.25ug/ml)
staining of
paraffin...

Target

SNX26 / TCGAP (N Terminus)

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Prostate shows cytoplasm staining of secretory cells in the glands. Recommended concentration, 1.25ug/ml.WB: Preliminary experiments gave an approx 35kDa band in human brain (hippocampus) lysates after 1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 121kDa according to NP_443180.2. The 35kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

ARSTDSLDPGE

Clonality

Polyclonal

MW

121

Entrez

115703

Dilution Range

ELISA: 1:32000, WB: 1 µg/ml, IHC-P: 1.25µg/ml

Expiration Date

12 months from date of receipt.