

Product Datasheet

SNAIL antibody (orb19327)

Description	Goat polyclonal antibody to SNAIL
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Snail homolog 1 / SNAIL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.1-0.3 µg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Preliminary experiments in lysates of cell lines HEK293, A431 and HepG2 gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	RKPSDPNRKPNY
Clonality	Polyclonal
MW	29.1
Dilution Range	ELISA: 1:164000, WB: 0.1-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse (A) and R...