

## Product Datasheet

### PABPC1L2B Antibody (C-term) (orb1932720)

|                            |                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Affinity Purified Rabbit Polyclonal Antibody (Pab)                                                                                                                             |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                         |
| <b>Reactivity</b>          | Human                                                                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                   |
| <b>Tested Applications</b> | WB                                                                                                                                                                             |
| <b>Target</b>              | This PABPC1L2B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 172-200 amino acids from the C-terminal region of human PABPC1L2B. |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.   |
| <b>Note</b>                | For research use only                                                                                                                                                          |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                     |
| <b>MW</b>                  | 22799 Da                                                                                                                                                                       |
| <b>Uniprot ID</b>          | <a href="#">Q5JQF8</a>                                                                                                                                                         |
| <b>NCBI</b>                | <a href="#">NP_001035971.1</a> , <a href="#">NP_001012995.1</a>                                                                                                                |
| <b>Dilution Range</b>      | WB: 1:1000                                                                                                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                |