


Product Datasheet

41883 antibody (orb19328)

Description	Goat polyclonal antibody to 41883
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Septin 1 / SEPT1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 45kDa band observed in lysates of cell line Jurkat (calculated MW of 42.0kDa according to NP_443070.1). Recommended concentration: 2-4 µg/ml.
Protein Sequence	GARDRASRSKLSRQ
Clonality	Polyclonal
MW	42.4
Entrez	1731
Dilution Range	ELISA: 1:32000, WB: 2-4 µg/ml
Expiration Date	12 months from date of receipt.



orb19328
(2ug/ml)
staining of
Jurkat lys...