

Product Datasheet

USP39 Antibody (Center) (orb1932936)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This USP39 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 237-265 amino acids from the Central region of human USP39.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65381 Da
Uniprot ID	Q53GS9
NCBI	NP_001243656.1, NP_001243654.1, NP_001243657.1, NP_006581.2, NP_001243655.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.