

Product Datasheet

SEPT6 Antibody (C-term) (orb1932937)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This SEPT6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 400-426 amino acids from the C-terminal region of human SEPT6.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49717 Da
Uniprot ID	Q14141
NCBI	NP_665799.1 , NP_665798.1 , NP_665801.1 , NP_055944.2
Dilution Range	WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.