

Product Datasheet

CORO1C Antibody (C-term) (orb1932939)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This CORO1C antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 377-403 amino acids from the C-terminal region of human CORO1C.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53249 Da
Uniprot ID	Q9ULV4
NCBI	NP_055140.1 , NP_001098707.1 , NP_001263400.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.