

Product Datasheet

ALS2CR8 Antibody (N-term) (orb1932940)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This ALS2CR8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 132-158 amino acids from the N-terminal region of human ALS2CR8.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	80698 Da
Uniprot ID	Q8N187
NCBI	NP_001098056.1, NP_079020.13
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.