

Product Datasheet

HHCM Antibody (C-term) (orb1933024)

Description Affinity Purified Rabbit Polyclonal Antibody (Pab)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Target This HHCM antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 341-368 amino acids from the C-terminal region of human HHCM.

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

Uniprot ID [Q05877](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.