

Product Datasheet

HIST1H2BH/HIST1H2BK/HIST3H2BB Antibody (Center) (orb1933033)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This HIST1H2BH/HIST1H2BK/HIST3H2BB antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 26-52 amino acids from the Central region of human HIST1H2BH/HIST1H2BK/HIST3H2BB.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	13892 Da
Uniprot ID	Q93079
NCBI	NP_003515.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.