

Product Datasheet

POMZP3 Antibody (N-term) (orb1933163)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This POMZP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human POMZP3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20620 Da
Uniprot ID	Q6PJE2
NCBI	NP_694537.1 , NP_036362.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.