

Product Datasheet

PCDHGA5 Antibody (N-term) (orb1933205)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This PCDHGA5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 254-280 amino acids from the N-terminal region of human PCDHGA5.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	100935 Da
Uniprot ID	Q9Y5G8
NCBI	NP_061741.1, NP_114443.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.