

Product Datasheet

NHEJ1 antibody (orb19337)

Description

Goat polyclonal antibody to NHEJ1

Species/Host

Goat

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Target

XRCC4-like factor / NHEJ1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 38kDa band observed in lysates of Human Skin, Testis and Thyroid gland (calculated size of 33.3kDa according to NP_079058.1). The observed molecular weight corresponds to findings with antibodies from other sources. The 38kDa band was successfully blocked by incubation with the immunizing peptide. Recommended concentration 0.1-0.3 µg/ml.

Protein Sequence

QRPQLSKVKRKKPR

Clonality

Polyclonal

MW

33.3

Entrez


79840

Dilution Range


ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 2.5-3.8µg/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of Human Thyroid L...



Immunohistochemical staining of Human Pl...