

Product Datasheet

ZMYM6 Antibody (N-term) (orb1933823)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This ZMYM6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 183-210 amino acids from the N-terminal region of human ZMYM6.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	148089 Da
Uniprot ID	O95789
NCBI	NP_009098.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.