

Product Datasheet

Mouse Rps6ka1 Antibody (C-term) (orb1934028)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Mouse Rps6ka1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 556-583 amino acids from the C-terminal region of mouse Rps6ka1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	81595 Da
Uniprot ID	P18653
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.