

Product Datasheet

IL20RA Antibody (Center) (orb1934029)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This IL20RA antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 141-169 amino acids from the Central region of human IL20RA.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62485 Da
Uniprot ID	Q9UHF4
NCBI	NP_001265651.1, NP_055247.3, NP_001265652.1, NP_001265653.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.