

Product Datasheet

RBM23 Antibody (Center) (orb1934123)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This RBM23 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 128-157 amino acids from the Central region of human RBM23.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48731 Da
Uniprot ID	Q86U06
NCBI	NP_001070820.1 , NP_060577.3 , NP_001070819.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.