

Product Datasheet

KLHL12 Antibody (C-term) (orb1934125)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This KLHL12 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 376-404 amino acids from the C-terminal region of human KLHL12.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	63277 Da
Uniprot ID	Q53G59
NCBI	NP_067646.1
Dilution Range	WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.