

Product Datasheet

Mouse Stk32a Antibody (C-term) (orb1934136)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Mouse Stk32a antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 347-374 amino acids from the C-terminal region of mouse Stk32a.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46509 Da
Uniprot ID	Q8BGW6
NCBI	NP_848864.1
Dilution Range	WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.